



Precise Numerical Results

Aflatoxin M₁, a toxin found in milk, is derived from ingested aflatoxin. A known carcinogen, the detection of M₁ is critical because this aflatoxin is found in milk, a major source of nutrition for infants, children, and adults. Throughout the world, there are significant regulations on the presence of M₁ in dairy products. VICAM's Afla M₁™ HPLC is a fast flow column that delivers recovery rates from 75% to 95% depending on the sample.

Using monoclonal antibody based affinity chromatography, Afla M₁ HPLC is an HPLC only test for the detection of aflatoxin M₁ using immunoaffinity columns. Performed in less than 30 minutes, Afla M₁ HPLC can measure levels of Afla M₁ as low as 10 ppt and as high as 3 ppb.

An ideal cleanup step for any HPLC, Afla M₁ HPLC has a long shelf life and provides precise numerical results in parts per billion or parts per trillion levels. The preferred test of laboratories around the world, no other test comes close for speed, quantification, and specificity.



HERE'S HOW SIMPLE AFLA M₁ HPLC IS TO USE:

PREPARE SAMPLE

- Centrifuge

ABSORB AND ELUTE

- Pass sample over Afla M₁ HPLC affinity column
- Wash column
- Elute aflatoxin M₁ and collect in cuvette

MEASURE

- Inject eluate into HPLC
- Determine aflatoxin

BENEFITS

Exclusive – Specifically for HPLC use

Durable – Long shelf life; requires no refrigeration

Versatile – Can be used with a variety of milk samples

Quick – Sample ready in less than 30 minutes

Wide Range – Detects levels as low as 10 ppt and as high as 3 ppb

Safe – Requires less toxic materials than other methods



ORDERING INFORMATION

Cat. No.	Description	Cat. No.	Description
G1007	Afla M ₁ ™ HPLC Columns (25 per box)	20700	Wash Bottle, 500 mL
G1106	Afla M ₁ HPLC Single Position Pump Stand, 10 mL & 60 mL*	21010	Cuvette Rack
G1108	Cuvette Holder	23040	Vortex Mixer
G1118	WB Column Coupling (6 per pack)	34000	Disposable Cuvettes (250 per pack)
G1130	Acetonitrile, HPLC Grade (4 x 4 L)	35016	Methanol, HPLC Grade (4 x 4 L)
G4033	Micro pipettor, 1 mL	36010	Disposable Plastic Beakers (25 per pack)
20050	Graduated Cylinder, 50 mL		
20100	Digital Scale with AC Adapter		
20250	Graduated Cylinder, 250 mL		

Also required; must be sourced locally:

Centrifuge
Distilled, Reverse Osmosis, or Deionized Water

* Multiple position pump stands available

Subject to change without notice.

© 2010 Waters Corporation. Waters, The Science of What's Possible, VICAM, AflaM₁ are trademarks of Waters Corporation.

The analytical methods presented in this data sheet have been researched and developed by Waters to be used exclusively with AflaM₁ products. These analytical methods have been validated in the Waters laboratories to perform to the specifications indicated in the AflaM₁ procedures. The user assumes all risk in using AflaM₁ analytical procedures and products. Waters makes no warranty of any kind, expressed or implied, other than that AflaM₁ products conform to Waters' printed specification and quality control standards. Waters will, at its option, repair or replace any product, or part thereof, which proves to be defective in workmanship or material. VICAM's undertaking to repair or service such products is exclusive and is in lieu of all other warranties whether written, oral, expressed, or implied, including any implied warranty of merchantability or fitness for a particular purpose. Waters shall have no liability for anticipated or lost profits or any loss, inconvenience or damage whether direct, indirect, incidental, consequential or otherwise, to person or property, or for strict liability or negligence arising from in connection with the use of these assay procedures or AflaM₁ products.

Orders:

1848 N Deffer Drive
Nixa, MO 65714 USA
Tel: +1.877.228.4244
+1.417.725.6588
Fax: +1.417.725.6102
orders@vicam.com

Headquarters:

34 Maple Street
Milford, MA 01757 USA
Tel: +1.800.338.4381
+1.508.482.4935
Fax: +1.508.482.4972
www.vicam.com